

## TruNeb Reusable Nebulizer Kit

HCPCS:



- The TruNeb Reusable Nebulizer Kit is high-quality and provides patients with what they need for effective aerosol therapy treatments.
- It includes an angled mouthpiece, 7' supply tube, 10cc medication cup, and a reusable nebulizer.
- Delivers fast and effective treatment efficiently
- Easy-seal, threaded cap and anti-spill jar prevent leakage of costly medication
- Provides 80% of particles in respirable range
- Nebulizes in upright positions of at angles up to 45 degrees
- Dishwasher safe
- Includes an angled mouthpiece, 7' supply tube, high capacity medication cup (10cc), and a reusable nebulizer
- Quantity: 10/case
- For a high-quality nebulizer kit, purchase the TruNeb Reusable Nebulizer Kit.
- Disposable option also available.

### PRODUCTS

### Description

50016

TRUNEB REUSABLE NEBULIZER KIT 10EA/CV