

AirLift Oxygen Tank Carrier, D/E Cylinders

HCPCS:



- The Airlift Oxygen Tank Carrier makes transporting your oxygen cylinder easy and comfortable.
- It includes a clear vinyl window to conveniently check cylinder's contents gauge without taking the cylinder out of the bag.
- Comfort flex shoulder strap and convenient hand grip keep oxygen tank transport easy and comfortable
- Can easily attach to wheelchair or scooter
- Made of non-flammable, durable, and non-abrasive materials; machine washable
- Includes dual zippers for ease of access and an accessory pocket
- Can hold D or E oxygen tanks
- Color: Navy
- To conveniently and comfortably transport your oxygen tank, purchase the AirLift Oxygen Tank Carrier.

PRODUCTS	Description
44N	AIRLIFT WHEELCHAIR CYLINDER BAG 50EA/CV