

Aluminum D Cylinder with Toggle Valve, 6/PK

HCPCS:

- The Catalina D Cylinder with CGA-870 Toggle Valve holds up to 425 liters of oxygen.
- This is available with wrench, toggle, or dual control valves in pallet or 6-pack quantities.
- High strength aluminum alloy up to 40% lighter than steel brushed body
- Green dome and heat sensitive clear coat finish.
- Easy to clean and maintain
- Includes toggle valve; also available with wrench or dual control valves
- Private labeling options available with your choice of stamping or silk screen.
- Pallet minimums are required, additional fees may apply.
- Individual weight: 6.32 lbs
- Quantity: 6/pack
- Purchase your Cylinders by the pallet and get free shipping on your cylinder order.

PRODUCTS

Description

PX-8702-1T	D Cylinder, toggle valve, 6/PK 6EA/PKG
------------	--
