

## Aluminum E Cylinder with Toggle Valve, 6/PK

HCPCS:

- We offer empty aluminum oxygen cylinders in standard sizes packaged to meet your requirements.
- The E Cylinder with CGA-870 Toggle Valve holds up to 679 liters of oxygen.
- This is available with wrench, toggle, or dual control valves in pallet or 6-pack quantities.
- High strength aluminum alloy up to 40% lighter than steel brushed body
- Green dome and heat sensitive clear coat finish.
- Easy to clean and maintain
- Includes toggle valve; also available with wrench or dual control valves
- Private labeling options available with your choice of stamping or silk screen.
- Pallet minimums are required, additional fees may apply.
- Individual weight: 8.5 lbs
- Quantity: 6/pack

### PRODUCTS

### Description

---

|            |  |
|------------|--|
| PX-8703-1T | E Cylinder, toggle valve, 6/PK 6EA/PKG |
|------------|--|

---