

Aluminum M60 Cylinder with Spigot Valve, **PALLET**

HCPCS:

- We offer empty aluminum oxygen cylinders in standard sizes packaged to meet your requirements.
- The Catalina M60 Cylinder with CGA-540 Spigot Valve holds up to 1699 liters of oxygen.
- This is available in pallet or 6-pack quantities.
- High strength aluminum alloy up to 40% lighter than steel brushed body
- Green dome and heat sensitive clear coat finish.
- Easy to clean and maintain
- Includes carrying handle
- Private labeling options available with your choice of stamping or silk screen.
- Pallet minimums are required, additional fees may apply.
- Pallet quantity: 42 cylinders per pallet
- Purchase your Cylinders by the pallet and get free shipping on your cylinder order

PRODUCTS	Description
PX-5405-1B	M60 Cylinder, spigot valve, 1EA/