

## Aluminum M9/C Cylinder with Wrench Valve, 6/PK

HCPCS:

- We offer empty aluminum oxygen cylinders in standard sizes packaged to meet your requirements.
- The Catalina M9/C Cylinder with CGA-870 Wrench Valve holds up to 246 liters of oxygen.
- This is available with wrench, toggle, or dual control valves in pallet or 6-pack quantities.
- High strength aluminum alloy up to 40% lighter than steel brushed body
- Green dome and heat sensitive clear coat finish.
- Easy to clean and maintain
- Includes wrench valve; also available with toggle or dual control valves
- Private labeling options available with your choice of stamping or silk screen.
- Pallet minimums are required, additional fees may apply.
- Individual weight: 4.65 lbs
- Quantity: 6/pack

### PRODUCTS

### Description

---

PX-8701-1W

M9/C Cylinder, wrench valve, 6EA/PKG

---