

Clear Plastic Nebulizer/Suction Pump Bag, 1 mil

HCPCS:



- This Clear Plastic Nebulizer/Suction Pump Bag is ideal for homecare and long-term care environments.
- It meets the accreditation requirements by various accreditation organizations such as JCAHO, ACHC, and CHAPS by keeping contaminated equipment separate from clean, ready-to-deliver equipment.
- 15" x 9" x 24", fits split spring beds
- Clear plastic material
- Thickness: 1 mil
- Quantity: 500/roll
- To best protect your patient's nebulizer or suction pump, purchase your Clear Plastic Nebulizer/Suction Pump Bag.
- Compass Health Brands offers a variety of equipment bags and covers for all your medical care needs.

PRODUCTS

Description

BAG-15924

NEB/SUCTION EQUIP BAGS 15X9X24