

Clear Plastic Patient Set-Up Bag, 16x12

HCPCS:



- This Clear Plastic Patient Set-Up Bag is ideal for homecare and long-term care environments.
- It meets the accreditation requirements by various accreditation organizations such as JCAHO, ACHC, and CHAPS by keeping contaminated equipment separate from clean, ready-to-deliver equipment.
- 16" x 12", keeps respiratory equipment clean and dry for transport and storage
- Clear plastic material with blue print and space for labeling and patient identification
- Includes 4" gusset drawstring bag
- Quantity: 500/box
- To best protect your patient's respiratory equipment, purchase your Clear Plastic Patient Set-Up Bag.
- Compass Health Brands offers a variety of equipment bags and covers for all your medical care needs.

PRODUCTS

Description

BAG-SETUP

PATIENT SET UP BAG, 16X12