

Disposable M6 Cylinder Tote, Holds 6, 10/Pack

HCPCS:

- Lightweight, versatile cardboard carrier simplifies patient set-up, delivery and storage of oxygen cylinders.
- For use with M6 size medical oxygen cylinders
- Heavy-duty cardboard cylinder tote
- Lightweight, moisture-resistant construction
- Easy grip handle for comfortable carry, ideal for patients transporting their own oxygen cylinders
- Assembles in a snap! Open, secure center walls, push in corner tabs and go!
- Designed for patient or office use
- Holds six (6) M6 cylinders
- Sold as 10 pack

PRODUCTS

Description

150-2005PK10

M6 CYLINDER TOTE - HOLDS 6 10EA/PKG