

Equipment Cover, 1.5 mil, 24x30

HCPCS:

- This Clear Plastic Equipment Cover is ideal for homecare and long-term care environments.
- It meets the accreditation requirements by various accreditation organizations such as JCAHO, ACHC, and CHAPS by keeping contaminated equipment separate from clean, ready-to-deliver equipment.
- 24" x 30" general equipment cover
- Clear plastic material
- Thickness: 1.5 mil
- Quantity: 500/roll
- Box size: 26" x 10" x 10"
- To best protect your homecare equipment, purchase your Clear Plastic Equipment Cover.
- Compass Health Brands offers a variety of equipment bags and covers for all your medical care needs.

PRODUCTS

Description

BAG-2430R

GEN PURPOSE EQUIP BAGS 1.5 MIL 24X30