

Equipment Cover, 1.5 mil, 26x36

HCPCS:

- This Clear Plastic Equipment Cover is ideal for homecare and long-term care environments.
- It meets the accreditation requirements by various accreditation organizations such as JCAHO, ACHC, and CHAPS by keeping contaminated equipment separate from clean, ready-to-deliver equipment.
- 26" x 36" general equipment cover
- Clear plastic material
- Thickness: 1.5 mil
- Quantity: 250/roll
- Box size: 27" x 8" x 8"

PRODUCTS

Description

BAG-2636R

GEN PURPOSE EQUIP BAGS 1.5 MIL 26X36