

Equipment Cover, 1.5 mil, 30x42

HCPCS:

- This Clear Plastic Equipment Cover is ideal for homecare and long-term care environments.
- It meets the accreditation requirements by various accreditation organizations such as JCAHO, ACHC, and CHAPS by keeping contaminated equipment separate from clean, ready-to-deliver equipment.
- 30" x 42" general equipment cover
- Clear plastic material
- Thickness: 1.5 mil
- Quantity: 250/roll
- Box size 32" x 8" x 8"

PRODUCTS	Description
BAG-3042R	GEN PURPOSE EQUIP BAGS 1.5 MIL 30X42