

Oxygen Tubing

HCPCS: E0431



- Kink resistant 3-channel extension oxygen tubing comes in a variety of lengths from 7 to 50 feet
- Compatible with most oxygen tubing rigid, swivel, and Y connectors, as well as inline concentrator water traps
- Latex free

PRODUCTS	Description
TUB NK 07	No-kink oxygen tubing 7 feet,50/cs
TUB NK 14	No-kink oxygen tubing 14 ft, 50/cs
TUB NK 25	No-Kink Oxygen Tubing,25 ft, 25/cs
TUB NK 40	No-kink oxygen tubing 40ft, 20/cs
TUB NK 50	No-kink oxygen tubing 50 feet, 20/cs