

iGo2 Portable Oxygen Concentrator with Bluetooth

HCPCS: E1390 + E1392




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

The iGo2 POC provides patients additional peace of mind and DME Providers the ability for Remote Fleet Monitoring.

- Settings 1-5 Pulse Dose
- SmartDose Auto-Adjusting Converter Technology included
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive converter trigger on the POC market
- 8.4" (H) x 3.5"(W) x 8.6"(L)
- Bluetooth Connectivity
- SmartLinkO2 App
- Remote Fleet Monitoring
- SmartLinkO2 Cloud
- Assembled in the USA
- 5 yr warranty, 3 yr compressor, 1 yr sieve beds & battery

PRODUCTS

Description

125D-BT	iGO2 Port O2 Conctrtr w BT
125D-BT-XB	iGO2 Port O2 Conctrtr w BT & X battry

SPECIFICATIONS

Description

CA Prop 65	no
------------	----

ACCESSORIES

Description

125D-613	iGO2 Battery pack, EA
125D-670	iGO2 Carrying case, EA
DV51D-606	US POWER CORD, EA
DV68-620	iGO2 AC power supply, EA
125CH-613	iGO2 Battery Charger (US), EA
DV6X-619	DC car adaptor for IntelliPAP 2, EA