

Roscoe 25' Green Supply Tubing, Kink Resistant

HCPCS:



- Roscoe's Oxygen Supply Tubing provides patients with high-quality performance and maximum safety.
- It is latex-free and available in clear or a distinctive green tint color for greater visual contrast and improved safety.
- Three-channel inner lumen design is kink- and crush-resistant, providing uninterrupted flow of oxygen for maximum safety
- Tubing lies out straight when uncoiled, minimizing "torque" or pull on cannula
- Universal connectors at each end come standard
- Latex free
- Length: 25 feet
- Green-tinted for greater visual contrast; also available in clear
- Quantity: 25/case
- Not made with natural rubber latex
- For high quality tubing for your oxygen therapy, purchase Roscoe's Oxygen Supply Tubing.

PRODUCTS

Description

TUB-ROS25G

OXYGEN TUBING 25' GREEN 25EA/CV