

Roscoe 7' Clear Supply Tubing, Kink Resistant

HCPCS:



- Roscoe's Oxygen Supply Tubing provides patients with high-quality performance and maximum safety.
- It is latex-free and available in clear or a distinctive green tint color for greater visual contrast and improved safety.
- Three-channel inner lumen design is kink- and crush-resistant, providing uninterrupted flow of oxygen for maximum safety
- Tubing lies out straight when uncoiled, minimizing "torque" or pull on cannula
- Universal connectors at each end come standard
- Latex free
- Length: 7 feet
- Clear; also available in green tint color
- Quantity: 50/case
- Not made with natural rubber latex
- For high quality tubing for your oxygen therapy, purchase Roscoe's Oxygen Supply Tubing.

PRODUCTS	Description
TUB-ROS7	OXYGEN TUBING 7' CLEAR 50EA/CV