

Roscoe Adult Curved Soft Nasal Cannula with 4' Tubing

HCPCS:



- The Roscoe Nasal Cannulas are comfortable, economical and latex-free.
- They are made of extra-soft material, allowing them to fit comfortably over the ears.
- Soft and crush resistant latex-free material provides comfort and maximum safety
- Supply tubing lies straight when uncoiled to minimize pull on the cannula
- Supply tubing has six-channel inner lumen design that provides uninterrupted flow of oxygen
- Tubing and cannula are disposable
- Tubing length: 4 feet
- Quantity: 50/case
- Disposable and intended for single patient use; should be replaced every 30 days for best performance
- Not made with natural rubber latex
- Purchase your Roscoe Adult Curved Nasal Cannula for maximum comfort and safety at a low cost.

PRODUCTS

Description

CAN-ROS4

OXYGEN SOFT CANNULA W/ 4' TUBE 50EA/CV