

Roscoe Empty Humidifier Jar, 6psi

HCPCS:



- The Roscoe Empty Humidifier Jar is a single patient use product designed for long-term use and reliable performance.
- It operates quietly and will activate an audible pressure safety alarm at 6 psi.
- Large, easy-to-grip wing nut connector assures a secure connection to oxygen supply source
- Easy-on/-off lid provides a positive seal, no gasket required
- Jar permanently etched with large print minimum/maximum water level markings that are easy to read and will not rub or wash off
- Soft plastic diffuser resists clogging due to mineral build up from water
- Ideal for use with all oxygen sources including concentrators, liquid, tanks, and wall outlets
- 350 cc bottle
- Disposable, single-use
- Quantity: 50/case
- Not made with natural rubber latex

PRODUCTS	Description
HUM-ROS6	BUBBLE HUMIDIFIER WITH 6 PSI 50EA/CV