

Salter 3-Channel Oxygen Supply Tubing, 25 Foot, Clear

HCPCS:

- The Salter 3-Channel Oxygen Supply Tubing provides patients with an uninterrupted flow of oxygen for utmost safety.
- It is kink-resistant and crush-proof with rounded and tapered edges to prevent snagging on other surfaces.
- 3-channel design is kink- and crush-resistant, allowing for a safe, uninterrupted flow of oxygen
- Made of soft, comfortable, and disposable material
- Ribbed body allows for easy gripping when making connections
- Lies straight when uncoiled to minimize pull on the cannula
- Rounded and tapered edges won't snag or catch on other surfaces
- 25' length, available in a variety of other lengths
- Clear
- Quantity: 25/case
- For safe and reliable oxygen therapy, purchase the Salter 3-Channel Oxygen Supply Tubing.

PRODUCTS

Description

2025