

Salter 3-Channel Oxygen Supply Tubing, 50 Foot, Green

HCPCS:

- The Salter 3-Channel Oxygen Supply Tubing provides patients with an uninterrupted flow of oxygen for utmost safety.
- It is kink-resistant and crush-proof with rounded and tapered edges to prevent snagging on other surfaces.
- 3-channel design is kink- and crush-resistant, allowing for a safe, uninterrupted flow of oxygen
- Made of soft, comfortable, and disposable material
- Lies straight when uncoiled to minimize pull on the cannula
- Ribbed body allows for easy gripping when making connections
- Rounded and tapered edges won't snag or catch on other surfaces
- 50" length, available in a variety of other lengths
- Green tinted for greater visual contrast
- Quantity: 20/case
- For safe and reliable oxygen therapy, purchase the Salter 3-Channel Oxygen Supply Tubing.

PRODUCTS	Description
2050G	SALTER 50' OXYGEN TUBING-GREEN 20EA/CV