

Salter Adult Demand Cannula, 7' Tubing

HCPCS:



- The Salter Adult Demand Cannula provides patients of dual port oxygen conservation delivery systems with comfort and convenience.
- It features a dual lumen and comes equipped with 3-channel, kink-resistant supply tubing for a safe and continuous supply of oxygen.
- Soft, curved prongs fit comfortably in nose
- Comfortable over-the-ear style allows patients freedom of movement without losing oxygen
- Dual lumen for use with dual port style oxygen conserver device
- 7', kink- and crush-resistant supply tubing provides safety and a continuous supply of oxygen
- Disposable and intended for single-patient use; should be replaced every 14 days for best performance
- Quantity: 25/case

PRODUCTS

Description

4807	SALTER DEMAND CANNULA W/7' TUBE 25EA/CV
------	---