

## Salter Adult Nasal Cannula with 4' Oxygen Supply Tubing

HCPCS:

- The Salter Adult Nasal Cannula provides patients with comfort and convenience.
- It comes equipped with 3-channel, kink-resistant supply tubing for a safe and continuous supply of oxygen.
- Soft, curved prongs fit comfortably in nose
- Comfortable over-the-ear style allows patients freedom of movement without losing oxygen
- 4', kink- and crush-resistant supply tubing provides safety and a continuous supply of oxygen
- Disposable and intended for single-patient use; should be replaced every 14 days for best performance
- Quantity: 50/case
- For a comfortable and safe cannula option, purchase the Salter Adult Nasal Cannula with Supply Tubing.

### PRODUCTS

### Description

|        |                                  |
|--------|----------------------------------|
| 1600-4 | SALTER CANNULA W/4' TUBE 50EA/CV |
|--------|----------------------------------|