

Salter Adult Nasal Cannula with E-Z Wrap Foam Cushions & 7' Supply Tubing

HCPCS:

- The Salter Adult Nasal Cannula provides patients with comfort and convenience.
- It comes equipped with 3-channel, kink-resistant supply tubing for a safe and continuous supply of oxygen and E-Z Wrap Foam Cushions for added comfort.
- Soft, curved prongs fit comfortably in nose
- Comfortable over-the-ear style allows patients freedom of movement without losing oxygen
- Pre-attached E-Z Wrap foam cushions add comfort and are fully adjustable along the cannula tubes to ensure proper positioning for each patient
- 7', kink- and crush-resistant supply tubing provides safety and a continuous supply of oxygen
- Disposable and intended for single-patient use; should be replaced every 14 days for best performance
- Quantity: 25/case
- For a comfortable and safe cannula option, purchase the Salter Adult Nasal Cannula with Supply Tubing and E-Z Wrap Foam Cushions.

PRODUCTS	Description
1600TLC	SALTER TLC CANNULA W/7' TUBE 25EA/CV