

Salter E-Z Wraps

HCPCS:

- The Salter E-Z Wraps add comfort to cannulas to reduce pressure on skin and prevent ear chafing.
- They easily wrap around and slide over cannula tubing to provide a soft cushion around the ear.
- Soft foam adds cushioning that provides comfort and prevents ear chafing
- Easily slide over the head or face tube of an oxygen cannula, molding to tube
- Closed-cell foam repels moisture to best protect skin
- Quantity: 50/case
- To prevent air chafing typically associated with regular cannula use, purchase the Salter E-Z Wraps.

PRODUCTS	Description
1016	Hartway,W/C,16"R/Arm,Fxd Foot,1/cs