

## SALTER HIGH FLOW HUMIDIFIER

HCPCS:

- The Salter High-Flow Bubble Humidifier is a disposable, single-patient use device for use in health care facilities and home care.
- The humidifier is used to provide additional humidity to oxygen therapy for flow rates at 6 to 15 liters per minute (LPM).
- The Bubble Humidifier has a 6 PSI safety pop-off valve, which sounds and releases excess pressure if the oxygen tubing becomes blocked.
- Accommodates patients with oxygen flow rates between 6-15 LPM
- Large, easy-to-grip wing nut connector assures a secure connection to oxygen supply source
- 360-degree down spout diffuser increases overall humidity and keeps operation quiet
- Jar permanently etched with large print minimum/maximum water level markings that are easy to read and will not rub or wash off
- 6 PSI
- For superior humidification for patients with high oxygen flow levels, purchase the Salter High-Flow Humidifier.

---

### PRODUCTS

### Description

---

7900-0-25