

Salter Humidifier Adapter 15", Male Thread Grip, 50/CS

HCPCS:



- The Salter 3-Channel Oxygen Tubing provides patients with an uninterrupted flow of oxygen for utmost safety.
- It is kink-resistant and crush-proof with rounded and tapered edges to prevent snagging on other surfaces.
- 3-channel design is kink- and crush-resistant, allowing for a safe, uninterrupted flow of oxygen
- Made of soft, comfortable, and disposable material
- Includes male thread grip and ribbed connector for easy connections
- Lies straight when uncoiled to minimize pull on the cannula
- Rounded and tapered edges won't snag or catch on other surfaces
- 15" length, available in a variety of other lengths
- Clear
- Quantity: 50/case
- For safe and reliable oxygen therapy, purchase the Salter 3-Channel Oxygen Tubing.

PRODUCTS

Description

SO-676

SALTER HUMIDIFIER ADAPTER-15" 50EA/CV