

Salter In-Line Water Trap with Screw Ends

HCPCS:

- The Salter In-Line Water Trap collects water and moisture that may form during oxygen therapy.
- It traps the condensation from the system, preventing water from entering the patient's nose.
- For use with in-line concentrators
- Simple, durable design that minimizes hazard of snagging on furniture
- Screw-on end caps
- Clear for visual inspection of water build-up
- Quantity: 25/case
- To give patients the benefit of humidification therapy without condensation entering their nose, purchase the Salter In-Line Water Trap.

PRODUCTS	Description
7001-25	SALTER IN-LINE WATER TRAP 1EA/