

## Salter Style O2 Cannula, Adult High Flow with 7', 3 Channel Tube

HCPCS:

- High Flow O2 Nasal Cannula features include:
- Delivers oxygen flows up to 15 LPM
- Special three channel headset tubing with enhanced reservoir face piece
- Unique, rectangular style head set tubing slide bolo
- Seamless, extra soft, reservoir face piece anatomically contoured to fit upper lip
- Bright green supply tubing from the cannula to oxygen source

| <b>PRODUCTS</b> | <b>Description</b>                       |
|-----------------|--|
| 1600HF-7-25     | SALTER HI FLOW CANNULA W/7' TUBE 25EA/CV |