

Salter Twist-and-Pull Adapter, Barbed

HCPCS:

- The Salter Twist-and-Pull Adapter makes connecting any two lengths of standard oxygen tubing convenient and simple.
- It has a twist-and-pull body for easy connecting and barbed ends that hold tubing together securely once connected.
- Plastic adapter connects any two lengths of oxygen tubing with standard 5mm to 7mm end connectors
- Concentric rings make twist-and-pull connection easy
- Features raised barbed ends for tight and secure connection
- Clear to visually check for connector patency
- Quantity: 50/case
- To easily and securely connect your oxygen tubing, purchase the Salter Twist-and-Pull Adapter.

PRODUCTS	Description
1215	