

## Zephyr C Cylinder Shoulder Bag

HCPCS:

- The Zephyr Oxygen Cylinder Shoulder Bag allows you to conveniently transport your oxygen tank and system.
- It includes a mesh zippered window to easily view cylinder's contents gauge without taking the cylinder out of the bag.
- Adjustable shoulder strap minimizes stress from cylinder's weight
- Mesh zippered window for contents gauge visibility
- Made of heavy-duty polyester for durability and easy maintenance
- Includes dual side cannula portholes and dual zippers for ease of access
- Can hold C oxygen tanks
- For an easy and comfortable way to transport your C oxygen tank, purchase the Zephyr C Cylinder Shoulder Bag.

| PRODUCTS | Description                          |
|----------|--------------------------------------|
| 33       | ZEPHYR SHOULDER BAG FOR M9/C 60EA/CV |