

Suction Therapy Accessories

HCPCS: A7000, A7002, A9999

- Disposable Container Filter Cartridge Assembly for 800ml. Qty 12 per Package

PRODUCTS	Description
4650D-609	PUSH-IN FEET (4-PK)
7305D-605	SEP-SALES COMPRESSOR ASS'Y
7305D-606	SEP SALES CARRYING CASE
7305D-607	REPLACEMENT GAUGE
7305D-608	BACTERIA FILTER (NON-STERILE) 12EA/PKG
7305D-609	SEP-SALES HEAD CAP ASS'Y
7305D-610	CIRCUIT BOARD, 7305D SERIES
7305D-612	REPLACEMENT TUBING 6-PK
7305D-614	ElbowConct,Filter, 3.375" tubing, 1ea/ea
7305D-615	CUP SEAL (DFT) SEPARATE SALE PKG
7305D-616	BOT COVER ASSY (DFT) SEP SALE PACK
7305D-617	TOP COVER SEPARATE SALE PACKAGE
7305D-618	LED LABEL -DFT SEPARATE SALE PACKAGE
7305D-619	ELBOW (DFT) 6-PK SS
7305D-623	CYLINDER W/VALVES
7305D-632	Int Fltr Disp Canstr PK48
7305D-633	Int Filter Disp Canister Comp Kit Assmd
7305D-633-25	Int Fltr Disp Cnstr Comp Kit Assmd Cs 25
7305D-635	800ML DISP CONTR FILTER Cart. Asml,12 P
7305P-605	LED LABEL 5 PK
7305P-610	CIRCUIT BOARD, 7305P SERIES
7305P-614	SEP SALES BATTERY ASSY,1/ea
7305P-620	BATTERY RETROFIT KIT
7305P-631	CONTINENTAL EUROPE CORD