

## Full Face ComfortFit Deluxe CPAP Mask

HCPCS: A7030 A7031 A7035-1,A7030 A7031 A7035-1



FIGURE A



FIGURE B

- The ultimate in comfort and performance for PAP therapy
- Maximizes compliance by using a soft silicone cushion with frame stabilizer that redistributes pressure evenly over a large surface area
- The 100 FD series of masks conforms comfortably to the user's face without the need for forehead pad, allowing the user the ability to read, watch TV and feel less confined during PAP therapy

PRODUCTS	Description
100FDL	MASK,FULLFACE,COM-FIT DLX,LG
100FDL-NH	FULLFACE MASK,COM-FIT DLX,LGW/OHDG, 1/EA
100FDM	MASK,FULLFACE,COM-FIT DLX,MED
100FDS	MASK,FULLFACE,COM-FIT DLX,SM
100FDS-NH	FULLFACE MASK,COM-FIT DLX,SMW/OHDG, 1/EA
ACCESSORIES	Description
B10344-001	CUSHION,100FDM,100FD+M, EA
B10343-001	CUSHION,100FDS,100FD+S, EA
B10345-001	CUSHION,100FDL,100FD+L, EA
100FDL-NH	FULLFACE MASK,COM-FIT DLX,LGW/OHDG, 1/EA, EA
100FDM-NH	FULLFACE MASK,COM-FIT DLX,MDW/OHDG, 1/EA, EA
C10200-018	HEADGEAR,100FD,100FD+,100ND+, EA
100FDS-NH	FULLFACE MASK,COM-FIT DLX,SMW/OHDG, 1/EA, EA