

V2 Full Face Mask

HCPCS: A7030+A7035



- Headgear uses new quick-release strap clips and a foam material for improved comfort
- 360-degree mask swivel for expanded patient freedom
- Detachable swivel allows user to disconnect
- Single-molded mask frame and cushion with integrated chin support
- Anatomically contoured design with ribbed support to improve seal and comfort

PRODUCTS	Description
113485	V2 Fl Fce Msk w/ HdGr - Lg
SPECIFICATIONS	Description
Warranty	90 days
ACCESSORIES	Description
669181	V2 Fl Fce Msk Fce pc - Med, EA
201493	V2 Fl Fce Msk HdGr - Med, EA
201492	V2 Fl Fce Msk HdGr - Sm, EA
201494	V2 Fl Fce Msk HdGr - Lg, EA
113571	V2 Fl Fce Msk Fce pc w/Swivel - Med, EA