

Sleepnet Ascend Full Face Mask w/ EZFit Headgear

HCPCS:



- Sleepnet Ascend Full Face Mask System with EZFit Headgear - includes small, medium and large cushion for the perfect fit.
- Sleepnet's Ascend CPAP Nasal Mask Series features advanced cushion technology with the replaceable AIR gel cushion.
- For added convenience, patients can replace the cushion without having to replace the whole mask.
- The patented AIR gel provides maximum comfort and minimum weight helping to reduce pressure points and skin irritation for better patient compliance.
- The flexible spacebar has a unique bendable design that provides the patient with superior stability without contacting the skin.

PRODUCTS

Description

50825

SLEEPNET ASCEND FULL FACE CPAP 30EA/CV