

## AIRgel Nasal Mask Replacement Cushion, Large

HCPCS:

- Sleepnet's Ascend CPAP Mask Series features advanced cushion technology with the replaceable AIRgel cushion.
- For added convenience, patients can replace the cushion without having to replace the whole mask.
- The patented AIRgel provides maximum comfort and minimum weight helping to reduce pressure points and skin irritation for better patient compliance.
- The flexible spacebar has a unique bendable design that provides the patient with superior stability without contacting the skin.
- Detachable cushion allows you to replace your cushion without replacing your entire mask.
- Flexible spacebar is designed to provide superior stability without contacting the skin
- Active headgear connector moves when the patient moves to help maintain an effective seal
- Available in three sizes (small, medium and large) with simple and breathable quick-release EZFit headgear
- 360-degree Rotating Swivel for optimal tube positioning and easier mobility
- Patented AIRgel with an ergonomic design and luxurious feel reduces pressure while maintaining an effective seal
- Available in your choice of Full Face or Nasal Mask

| PRODUCTS | Description                     |
|----------|---------------------------------|
| 50203    | ASCEND NASAL CUSHION LG 50EA/CV |