

DreamEasy Nasal Mask with Headgear, Medium

HCPCS:

- The Roscoe Medical DreamEasy? Nasal CPAP Mask with Comfort Cushion? is extremely lightweight and comes complete with headgear.
- The Comfort Cushion? seals comfortably around the nose and facial contours to improve overall compliance and minimize air leaks.
- The headgear buckle comes with a quick snap design making it easy to fasten and unfasten.
- Size: Medium
- Flexible, lightweight design allows for freedom of movement while maintaining a secure fit.
- Headgear with Velcro fasteners provide maximum versatility and movement.
- 360° rotational elbow with standard 22 mm tubing connector.
- The Comfort Cushion? made with medical liquid silicone is completely odorless.
- Comfort Cushion? can be removed for easy cleaning.
- Ergonomic design is airtight.
- Vent hole provides effective CO2 washout.

PRODUCTS	Description
CPM-DENM	DREAMEASY NASAL MASK W/ HGEAR-MD 50EA/CV