

## Full Face ComfortFit Deluxe CPAP Mask

HCPCS: A7030 A7031 A7035-1,A7030 A7031 A7035-1



FIGURE A



FIGURE B

- The ultimate in comfort and performance for PAP therapy
- Maximizes compliance by using a soft silicone cushion with frame stabilizer that redistributes pressure evenly over a large surface area
- The 100 FD series of masks conforms comfortably to the user's face without the need for forehead pad, allowing the user the ability to read, watch TV and feel less confined during PAP therapy

| PRODUCTS    | Description                                  |
|-------------|----------------------------------------------|
| 100FDL      | MASK,FULLFACE,COM-FIT DLX,LG                 |
| 100FDL-NH   | FULLFACE MASK,COM-FIT DLX,LGW/OHDG, 1/EA     |
| 100FDM      | MASK,FULLFACE,COM-FIT DLX,MED                |
| 100FDS      | MASK,FULLFACE,COM-FIT DLX,SM                 |
| 100FDS-NH   | FULLFACE MASK,COM-FIT DLX,SMW/OHDG, 1/EA     |
| ACCESSORIES | Description                                  |
| B10344-001  | CUSHION,100FDM,100FD+M, EA                   |
| B10343-001  | CUSHION,100FDS,100FD+S, EA                   |
| B10345-001  | CUSHION,100FDL,100FD+L, EA                   |
| 100FDL-NH   | FULLFACE MASK,COM-FIT DLX,LGW/OHDG, 1/EA, EA |
| 100FDM-NH   | FULLFACE MASK,COM-FIT DLX,MDW/OHDG, 1/EA, EA |
| C10200-018  | HEADGEAR,100FD,100FD+,100ND+, EA             |
| 100FDS-NH   | FULLFACE MASK,COM-FIT DLX,SMW/OHDG, 1/EA, EA |