

Sleepnet Ascend Nasal Mask System with EZFit Headgear

HCPCS:



- Sleepnet Ascend Nasal Mask System with EZFit Headgear - includes small, medium and large cushion for the perfect fit.
- Sleepnet’s Ascend CPAP Nasal Mask Series features advanced cushion technology with the replaceable AIR gel cushion.
- For added convenience, patients can replace the cushion without having to replace the whole mask.
- The patented AIR gel provides maximum comfort and minimum weight helping to reduce pressure points and skin irritation for better patient compliance.
- The flexible spacebar has a unique bendable design that provides the patient with superior stability without contacting the skin.

PRODUCTS	Description
50174	ASCEND NASAL MASK W/ HEADGEAR 50EA/CV