

Silicone Tubing, 4", for Sunrise 7305D-023

HCPCS:

- This suction tubing is high-quality and durable, made of medical grade unreinforced silicone.
- It is compatible with the Sunrise 7305D-023 Suction Machines.
- Durable and flexible platinum cured medical grade unreinforced silicone tubing
- Latex free, BPA free, and animal derived ingredient free
- Length: 4"
- Quantity: 12/pack
- Compatible with Sunrise 7305D-023 Suction Machines

PRODUCTS	Description
TK-SUC4	Silicone Tubing for Suction 4 1PKG/