

Suction Setup Kit

HCPCS:

- The Suction Setup Kit allows users to conveniently replace all necessary parts in their suction machine at once.
- It includes a yankauer with bulb tip and control vent, replacement elbow, bacteria filter, 800ml reusable collection jar, 6' tubing and 14" vinyl tubing.
- Each item in the kit can also be purchased separately if desired:
 - Bulb-tip Yankauer - #55K86V
 - Elbow - #ROS-610E
 - Bacteria Filter - #ILB-38
 - Collection Jar - #1160B-1
 - 6' Tubing - #55N56A
 - 14" Tubing - #196472-14
- To keep your portable suction machine performing optimally, replace its parts regularly and easily with the Roscoe Suction Kit.

PRODUCTS	Description
KIT-019	SUCTION SETUP KIT 1EA/