

## Tri-Flo Suction Catheter, 14 French Looped

HCPCS:

- The Cardinal Health Tri-Flo Single Suction Catheters are designed to provide high-quality suctioning that is less traumatic to patients.
- The triangular eyes and beveled tip make the suctioning process more comfortable without sacrificing high suction efficiency.
- 14 French looped catheter; available in a variety of sizes and types
- Provides high suction efficiency
- Beveled tip and triangular eyes make suctioning more comfortable and less traumatic to patients
- Includes control port
- Sterile and individually wrapped for added safety; intended for single patient use
- Quantity: 50/case

PRODUCTS	Description
T60C	Tri-Flo Suction Catheter, 1CV/