

YANKAUER SUCTION HANDLE 50/CS

HCPCS:

- The Bulb-Tip Yankauer allows for precise suctioning without damaging the surrounding tissue.
- It features a slip-resistant handle and smooth inner surface for high-quality performance.
- Provides precise and effective suctioning during medical procedures
- Does not include control vent; also available with vent (55K86V)
- Features slip-resistant handles, smooth inner surface and uniform tip-to-base internal diameter for fast evacuation
- Clear, transparent construction for easy visualization of fluids
- Ribbed five-in-one connector; fits various sizes of suction tubing
- Sterile, for single patient use
- Quantity: 50/case

PRODUCTS	Description
55K86	YANKAUER SUCTION HANDLE W/O VENT 50EA/CV