

ROSCOE ASPIRATOR, HIGH VACUUM

HCPCS:



- The Roscoe Aspirator is the perfect solution for all your suction needs.
- It is simple to use and has a durable design to withstand the rigors of homecare use.
- It produces a maximum vacuum of 675 mmHg and features a high-quality vacuum regulator with anti-vibration setting, an anti-bacteria filter design helps prevent contamination, oil-less, permanent lubrication, and a collection bottle with overflow protection.
- Designed for easy portability at just 7.7 lbs, it's easy to carry and easy to clean.

PRODUCTS

Description

ASP-ROS

ROSCOE ASPIRATOR, HIGH VACUUM 1EA/