

## Premium Guard Gel Overlay

HCPCS:



FIGURE A

- Best in class gel overlay optimizes pressure redistribution and is now available in up to 500 lbs weight capacity
- 15 gel bladders cover entire patient contact area for the prevention and treatment of pressure ulcers
- Full length "head-to-toe" gel bladders minimize "bottoming out" and "pitting"
- New convoluted foam topper vastly improves patient comfort and air circulation
- 4 elastic corner straps secure overlay to mattress
- Urethane coated nylon cover is water resistant and vapor permeable
- Vinyl base is stain resistant and waterproof
- Bladders contain an aqueous non-toxic gel and polystyrene beads to minimize gel migration
- Meets flammability standard TB117

PRODUCTS	Description
14893	PremGuardGel/FoamOverlay34x76x3.5in,1/cs
14901	Prem.GuardGel/FoamOverlay42x76x3.5",1/cs
SPECIFICATIONS	Description
Warranty	Foam: 18 Months Cover: 18 Months
Weight Capacity	300-500 lbs
HCPCS	E0185
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]
ACCESSORIES	Description
14893C	Prem.Guard Gel,3"Foam Overlay,Cover,1/ea, EA