

## Innerspring Mattress

HCPCS: E0271,E0271



- High quality innerspring design
- The unit has a 1" high-density, firm, combustion-modified fiber topper on each side
- Multi-layered with two 1" high-density, firm, combustion modified fiber toppers on the top and one 1" high density, firm, combustion-modified fiber topper on the bottom (15006EF only)
- Vinyl cover is waterproof and anti-static for easy use and care

PRODUCTS	Description
15006	Mattress 80" Innerspring, 1/EA
15006EF	Mattress,Inner Spring,ExtraFir,1/ea
15014	Mattress,Inner Spring36"x84"XF,1/ea

  

SPECIFICATIONS	Description
Warranty	1 Year
Weight Capacity	350 lb
Dimensions	80" x 36" x 6"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <a href="http://www.P65Warnings.ca.gov">www.P65Warnings.ca.gov</a> .]