

Innerspring Mattress

HCPCS: E0271,E0271



- High quality innerspring design
- The unit has a 1" high-density, firm, combustion-modified fiber topper on each side
- Multi-layered with two 1" high-density, firm, combustion modified fiber toppers on the top and one 1" high density, firm, combustion-modified fiber topper on the bottom (15006EF only)
- Vinyl cover is waterproof and anti-static for easy use and care

PRODUCTS	Description
15006	Mattress 80" Innerspring, 1/EA
15006EF	Mattress,Inner Spring,ExtraFir,1/ea
15014	Mattress,Inner Spring36"x84"XF,1/ea
SPECIFICATIONS	Description
Warranty	1 Year
Weight Capacity	350 lb
Dimensions	80" x 36" x 6"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]