

3" Foam Retail Cushion

HCPCS:



FIGURE A



FIGURE B




FIGURE C

- 3" foam cushion in attractive retail package with cutout to feel cover
- Comes with durable, machine washable, jacquard zippered cover
- Weight Capacity 250 lbs
- Made of soft, premium quality foam
- Versatile use for any location where seated comfort is needed

PRODUCTS

Description

RTL14910 Foam Cushion 18 x 16 x 3,1EA,12/CS

SPECIFICATIONS

Description

Warranty 1 Year Limited

Weight Capacity 250 lbs

Cover Material Jacquard Knit Cover