

Adult Anterior Safety Walkers

HCPCS: E1399-1,E1399-1



FIGURE A



FIGURE B



FIGURE C

- Height adjustable in 1" increments
- Adjustable-brake spring tension
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture
- Skid-resistant, non-marring, easy-rolling neoprene wheels
- Ideal for indoor and outdoor use
- Folds easily and stands on its own in the folded position (Figure A)

PRODUCTS

Description

CE 1000 BK	Adult SRoller,Silver Vein,1/cs
------------	--------------------------------

SPECIFICATIONS

Description

Warranty	Limited Lifetime
----------	------------------

Weight Capacity	400 lbs
-----------------	---------

Carton Shipping Weight	21 lbs.17 lbs.
------------------------	----------------

ACCESSORIES

Description

CE 1250	OxygenTank Holder,1/bx, BX
---------	----------------------------

CE 1500	Caster Lock Bracket,1pr/pkg, EA
---------	---------------------------------

CE 1400	Dial-A-Speed Tab,1pr/bx, BX
---------	-----------------------------

CE 1035 FP	ForearmPlatform&MountingBracket,1pr/bx, PR
------------	--

CE 1070L	Seat Harness,SRoller, Lg, 1/ea, EA
----------	------------------------------------