

## Support Kit

HCPCS: E1399-1,E1399-1



- Support Kit includes Trunk Harness, Abductor and Lateral Supports
- Trunk Harness:
  - Adjustable zippered trunk harness provides focused support on the sternum and rib cage
  - It's semi-stretch material allows controlled movement, prevents forward slumping and also maintains correct midline seated position
  - Straps with cam lock buckles secure the trunk harness
- Abductor:
  - Depth adjustable
  - Positions user's knee controls spastic muscles that might cause their legs to cross
- Lateral Supports:
  - Assists proper stabilization and alignment of the upper body
  - Height and width adjustable
  - Sold as a pair and each lateral is individually adjustable

### PRODUCTS

### Description

PRODUCTS	Description
FC 8025 LN	Support Kit, Lg, FC4000N, 1set/bx
FC 8025N	Support Kit, Sm, FC2000N, 1set/bx